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## Liebig-Dispenser

**DE:** Seite 2

---

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## Diffuseur Liebig

**FR:** page 4

---

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## Diffusore Liebig

**IT:** lato 6

---

Instructions of use for

## Liebig-Dispenser

**EN:** page 8

---

# Liebig-Dispenser

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
# Liebig-Dispenser

Effective and bee-safe treatment for Varroosis with formic acid 60% and/or 85%.

## Components

- 1 Base plate with reservoir and fixing prongs for a safe fixation.
- 2 Perforated paper-wick for easy and exact delivery.
- 3 Dropper with calibrated opening for continuous release.
- 4 Low profile bottle with graduation marks for easy filling and fast checking during treatments.

Designed for the use with the veterinarian medicines FORMIVAR® 60 and FORMIVAR® 85 for the treatment of Varroosis in honey bees (*Apis mellifera*).

## Mode of Action

Formic acid is delivered slowly and steadily through the calibrated opening in the dropper to the reservoir on the base plate. By capillary action the acid is distributed evenly over the surface of the paper-wick. The intensity of evaporation is controlled by the size of the paper-wick.

The paper is being folded along the fine perforated lines and excess paper is being removed. As long as there is formic acid in the bottle, evaporation will be constant. The treatment is finished as soon as the paper dries.

## Sizing of Paper-Wick

Respect the instructions printed on the paper-wick, green for formic acid 60% and red for formic acid 85%.

**Initial Sizing:** Determine hive type and size and the expected maximum temperatures for the first days of treatment. Consult your local weather forecast! At the point of intersection you see the recommended volume for the first or the second treatment. Reduce the size of a paper-wick to the recommended size by folding the paper along the perforated lines and by removing excess paper. Average daily evaporation can be determined dividing dosage by duration of treatment, e.g.

FORMIVAR® 60% for the first treatment 15–30 ml per day depending on hive size and 15–20 ml per day for the second treatment.

**Adjusting when evaporation is too high:** Is the daily consumption too high or you find the bees bearding up in front of the hive reduce the paper-wick(s) by removing one strip.

**Adjusting when evaporation is too low:** Add a second paper-wick of identical size on top of the already wet paper-wick increasing capillary effect. The paper-wick will soak faster and more formic acid will be transported to the edges.

## Instructions for a Treatment

### Prepare

### Dose

### Treat

1. Remove dropper.
2. Fill bottle to the recommended mark.
3. Insert dropper. Secure the bottle for transport and storage with child proof cap. Store at ambient temperature.
4. Size paper-wick respecting your hive size and type and expected maximum daily temperatures.
5. Place base plate in the center on top the brood frames and add paper-wick. Place the base plate horizontal and firm to prevent spilling of formic acid from reservoir.
6. Place bottle carefully without pressure with opening on top of the fixing prongs. Formic acid drips to the reservoir and is absorbed by the paper. **IMPORTANT:** Do not use chilled acid to avoid spilling due to acid warming to hive temperature.
7. Check flow of acid by slight tilting of bottle. Close the hive.
8. Close bottom board. Open flight entrance in commercial colonies. For nuclei leave flight opening unchanged to prevent robbery.
9. Check evaporated volume and adjust paper-wick if required at the earliest after one day. Once the bottle is empty and the paper dry remove Liebig-Dispenser from hive.

### Important:

- Keep dispenser horizontal.
- Do not use chilled formic acid to avoid spilling of acid due to acid warming to hive temperature.
- With expected high daily temperatures start treatment preferably in the early hours.
- To reduce stress for the bees do not start the treatment shortly before thunderstorms or heavy rains.
- Use only after honey harvest for best honey quality.
- Do not treat and feed simultaneously to prevent absorption of formic acid by fresh food and to prevent robbery.
- 1<sup>st</sup> treatment: End July to mid August before feeding for 3–7 days
- 2<sup>nd</sup> treatment: after completed feeding in September for 6–10 days
- Standard sizes for paper-wicks can be found on **page 10–11**. For larger volume hives we recommend to use two layers of paper-wicks to improve wetting of the paper with formic acid

## Cleaning/Additional Information

Remove propolis with HalaApi 898. The base plate may be cleaned from wax using APISODA in warm water. Do not clean the plate in the dishwasher. The base plate will bend at temperatures over 40 °C.

The treatment with formic acid and the Liebig-Dispenser is proven bee-safe for one and two storey hives and the Swiss hive. Open brood and emerging bees are very sensitive to formic acid.

A treatment with formic acid will always cause limited losses in brood. These losses have no negative impact on colony development. Queen losses will be avoided by following the instructions of use correctly.

Grundeinstellungen der Dochtpapierfläche für die Behandlung mit FORMIVAR® 60  
 Réglage du papier buvard pour un traitement avec FORMIVAR® 60  
 Impostazione di base della carta assorbente per il trattamento con FORMIVAR® 60  
 Sizing of paper-wick for a treatment with FORMIVAR® 60

FORMIVAR® 60%	CH-Kasten Ruche Suisse (Bürki) Arnica Svizzera	Nuclei Mini Plus Warré	Zander, DNM Langstroth Segeberger National Simplex, Smith 2×Warré	Dodant 2×Simplex AZ (Alberti- Žnidaršič) Commercial Deep National	2×Zander 2×DNM 2×Segeberger 2×Langstroth 4×Warré
	Vol.: 40 l	Vol.: ~20 l	Vol.: ~40 l	Vol.: ~50 l	Vol.: ~80 l
25–30 °C					
	1. 100 ml/7 d 2. 150 ml/10 d		1. 150 ml/7 d 2. 150 ml/10 d	1. 200 ml/7 d 2. 200 ml/10 d	
20–25 °C					
	1. 100 ml/7 d 2. 150 ml/10 d		1. 150 ml/7 d 2. 150 ml/10 d	1. 200 ml/7 d 2. 200 ml/10 d	
15–20 °C					
	1. 100 ml/7 d 2. 150 ml/10 d		1. 150 ml/7 d 2. 150 ml/10 d	1. 200 ml/7 d 2. 200 ml/10 d	

d = Tag/day/ Jour/giorno  
 1. = Erste Behandlung/first treatment/première traitement/primo trattamento  
 2. = Zweite Behandlung/second treatment/deuxième traitement/secondo trattamento  
 Vol. = Beutevolumen/Hive Volume/Volume de la ruche/volume d'arnia

Verwendung von FORMIVAR® 85% wird empfohlen.  
 Use of FORMIVAR® 85% is recommended.  
 L'utilisation de FORMIVAR® 85% est recommandée.  
 L'uso del FORMIVAR® 85% è consigliato.

Verwenden Sie 2 Dochtpapiere pro Behandlung.  
 Use 2 paper-wicks per treatment.  
 Utilisez 2 papiers par traitement.  
 Utilizzare 2 carte per trattamento.

Grundeinstellungen der Dochtpapierfläche für die Behandlung mit FORMIVAR® 85  
 Réglage du papier buvard pour un traitement avec FORMIVAR® 85  
 Impostazione di base della carta assorbente per il trattamento con FORMIVAR® 85  
 Sizing of paper-wick for a treatment with FORMIVAR® 85

FORMIVAR® 85%	CH-Kasten Ruche Suisse (Bürki) Arnica Svizzera	Nuclei Mini Plus Warré	Zander, DNM Langstroth Segeberger National Simplex, Smith 2×Warré	Dodant 2×Simplex AZ (Alberti- Žnidaršič) Commercial Deep National	2×Zander 2×DNM 2×Segeberger 2×Langstroth 4×Warré
	Vol.: 40 l	Vol.: ~20 l	Vol.: ~40 l	Vol.: ~50 l	Vol.: ~80 l
25–30 °C					
	1. 50 ml/7 d 2. 100 ml/10 d		1. 75 ml/7 d 2. 100 ml/10 d	1. 100 ml/7 d 2. 200 ml/10 d	
20–25 °C					
	1. 50 ml/7 d 2. 100 ml/10 d		1. 75 ml/7 d 2. 100 ml/10 d	1. 100 ml/7 d 2. 200 ml/10 d	
15–20 °C					
	1. 50 ml/7 d 2. 100 ml/10 d		1. 75 ml/7 d 2. 100 ml/10 d	1. 100 ml/7 d 2. 200 ml/10 d	

d = Tag/day/ Jour/giorno  
 1. = Erste Behandlung/first treatment/première traitement/primo trattamento  
 2. = Zweite Behandlung/second treatment/deuxième traitement/secondo trattamento  
 Vol. = Beutevolumen/Hive Volume/Volume de la ruche/volume d'arnia

Verwendung von FORMIVAR® 60% wird empfohlen.  
 Use of FORMIVAR® 60% is recommended.  
 L'utilisation de FORMIVAR® 60% est recommandée.  
 L'uso del FORMIVAR® 60% è consigliato.

Verwenden Sie 2 Dochtpapiere pro Behandlung.  
 Use 2 paper-wicks per treatment.  
 Utilisez 2 papiers par traitement.  
 Utilizzare 2 carte per trattamento.

Temperaturangaben entsprechen den zu erwartenden maximalen Temperaturen gemäß Wetterprognosen./Températures maximales attendues selon les prévisions météorologiques./La temperatura indicata è conforme alla temperatura massima prevista secondo le previsioni meteo./Temperature levels corresponds to the daily maximum temperatures according to weather forecast.

Temperaturangaben entsprechen den zu erwartenden maximalen Temperaturen gemäß Wetterprognosen./Températures maximales attendues selon les prévisions météorologiques./La temperatura indicata è conforme alla temperatura massima prevista secondo le previsioni meteo./Temperature levels corresponds to the daily maximum temperatures according to weather forecast.

1  
Tropfeinsatz entfernen  
Retirez le goutte-à-goutte  
Togliere il contagocce  
Remove dropper

2  
Ameisensäure dosieren  
Dosage de l'acide formique  
Dosare l'acido formico  
Measure formic acid

3  
Tropfeinsatz aufsetzen  
Remettre le goutte-à-goutte  
Riporre il contagocce  
Insert dropper

4  
Dochtpapier anpassen  
Réglez le papier buvard  
Adattare la carta assorbente  
Adjust size of paper-wick

5  
Dochtpapier auflegen  
Déposez le papier buvard  
Poggiare la carta assorbente  
Place paper-wick

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